




Accessories



Wall Mount Box, with Alarm Package

- Part Number: 5950
- To help stop theft, vandalism or misuse
- Box Features:
 - Size: 6 3/4 x 20 1/2 x 6
 - High impact crystal polystyrene of .110 wall thickness
 - Virgin acrylic Plaskolite® panel of .080 thickness
 - Chromium plated brass safety lock and key with weatherproof protective lock cover
 - Polypropylene plastic injected molded hammer with flexible torsion, ribbed grip surface, and clear nylon-jacketed stainless steel mechanical cable
- Alarm Features:
 - Highly visible “stop sign” design discourages wrongful use
 - Sounds a powerful 105 dB warning horn when activated
 - Deactivation pin provided for service by authorized personnel.
 - 9 VDC alkaline battery is included
 - Compact design 3.37 in. (83mm) in diameter and 1 in. (25mm) deep.
 - Three minute automatic reset upon activation.

Wall Mount Box Only

- Part Number: 5951

Alarm Only

- Part Number: 5952

Wall Mount Box Front Replacement

Panel Only

- Part Number: 5953
- Virgin acrylic Plaskolite® panel of .080 thickness

Mask Assembly

- Part Number: 5945
- Adult/Pediatric comfort fit mask
- 7' Kink-free, crush resistant hose

Wall Sign

- Part Number: B20559
- Easy to see, ‘V’ shaped projection-style sign
- White sign, Red Emergency Oxygen letters
- Size: 6” X 5”
- Custom signs are available to fit any décor

Display Unit

- Display unit available at a discounted rate
- See your sales representative for details